

BTSeq™

An NGS-based Innovative Sequencing Platform with High-fidelity and Low-cost

KEY BENEFITS

- **High fidelity sequencing**
 - High sequencing quality (NGS-based)
 - Specific and clear digitized sequencing result
- **Sequencing with longer reads**
 - Up to 20 kb in a single reaction
 - Short turnaround time (TAT), within a day.
- **No need for sequencing primers**
 - Primer information is not required.
 - Only the template is required for sequencing.
- **Low-cost sequencing**
 - Cost-effective compared with Sanger sequencing

The Workflow: Fast, Easy, High-fidelity and Low-cost Sequencing

In BTSeq™, DNAs are simultaneously fragmented and tagged with molecular barcodes by a simple enzyme reaction (Figure 1), which significantly reduces the experiment time. No extra equipment is necessary because of DNA fragments are made chemically and its simple workflow also decreases the experimental cost.

No limitation on DNA Samples

In BTSeq™, the DNA samples to determine their sequence are not limited to their type or length. Both PCR products, plasmids, and even the genomes of bacteria and viruses can be analyzed. Furthermore, sequencing unpurified or small amounts (~0.1 ng) of samples is available with high-quality. By covering the length from 150 bp up to 20 kb, BTSeq™ can analyze most DNAs that can be sequenced by Sanger method. Moreover, sequencing longer DNA is available in a single reaction. Therefore any consecutive (or repetitive) sequencing such as primer walking is not necessary. Since BTSeq™ doesn't require sequence-specific primers, it is possible to determine the sequence even if you don't have any available sequencing primers.


Introduction

BTSeq™ (Barcode-Tagged Sequencing) is a new sequencing platform developed to replace the existing methods, analyzing the DNA sequence rapidly and accurately without the DNA length limitation. It enables to prepare the sequencing-ready library for various length and types of DNA within 2 hours and to determine the DNA sequence using NGS and bioinformatics pipeline within a day. BTSeq™ is compatible with most existing NGS platforms, which is available for various applications and also can be automated easily for the high-throughput.

Accurate Calls

By using our proprietary bioinformatics pipeline after NGS, the sequencing reads are sorted according to the molecular barcodes, piled up to correct NGS errors, and finally completed the entire sequence of DNA. Therefore, BTSeq™ provides very accurate results with no ambiguity and does not need any reanalysis.

Figure 1. BTSeq™ Workflow



Developed a new enzyme and molecular barcodes, which simplifies sample preparation workflow and reduces turnaround time.

No limitation of samples
 - Unpurified sample
 - PCR product or plasmid
 - Sample amount (>0.1 ng)

Developed a bioinformatics pipeline to reduce NGS errors and to assemble longer consensus sequences up to 20 kbp.


Table 1. DNA sample information in the comparative study

Sample information
<ul style="list-style-type: none"> DNA type and length PCR product, n=801, 147 bp ~ 2.5 kb Plasmid, n=454, 2.7 kb ~12.9 kb
<ul style="list-style-type: none"> More than 80% of PCR product samples were not purified.
<ul style="list-style-type: none"> The range of sample concentration: 0.1ng/ul ~ 200 ng/ul 1 ul of samples (0.1 ng~200 ng) were used for BTSeq™.

Table 2. Summary of the comparison study of BTSeq™ and Sanger

# of samples	PCR product	Plasmid	Causes for discrepancy
DNA sample	801	454	
# of match	554	369	
# of discrepancy	247	85	
Causes from BTSeq™	36	0	Difficulty in the exact determination of long repeated sequences
Causes from Sanger	211	85	Poor sequencing results

Figure 3. The discrepancy of a Few Base Pairs due to Poor Quality Base Calls



Comparison of BTSeq™ and Sanger Sequencing

To validate the performance of BTSeq™ as an innovative sequencing platform, a variety of DNAs were sequenced using BTSeq™ and compared with their Sanger sequencing results. A total of 1,285 DNAs were selected, with varying types, lengths, quantities, and purity (Table 1). The whole and a part of respective DNAs were sequenced previously by the Sanger method and we compared both results. Since the starting parts and ending parts of Sanger data tend to be unclear, we only used the middle parts of them excluded from 100 bp region at starting parts and ending parts. After performing BTSeq™ using 1 ul from each DNA, the libraries were pooled and sequenced with paired-end 2×150 reads on the MiSeq system (Illumina). The NGS raw data were analyzed and assembled using bioinformatics pipeline of BTSeq™. Determined sequences were compared with those from Sanger sequencing.

Table 2 shows the result of BTSeq™. All of 1,285 DNA samples were sequenced clearly. Our bioinformatics pipeline produced the high-quality score for every sequence. Only a small amount of DNA (≥0.1 ng) is enough to analyze the various length from 150 bp to 13kb, regardless of its purity. When comparing the determined sequences with Sanger sequencing results (Table 2), 74% of the DNA sequences exactly matched to corresponding Sanger sequencing results. For 26% of DNAs that do not match each other, we analyzed in detail the Sanger chromatogram and the BTSeq™ quality score around the regions of the discrepancy. Poor Sanger sequencing results caused most discrepancies.

The most common pattern of discrepancies was when a few base pairs were mismatched. As shown in Figure 3, there was miscall due to the low-quality result in the Sanger chromatogram (ABI file), while the quality score of BTSeq™ base call for that sequence was high.

User-Friendly Reports

BTSeq™ provides the analysis report composed of a FASTA file and a chromatogram for the DNA sequence (Figure 2). It is similar to that of Sanger sequencing, so very familiar to most customers.

Figure 2. A BTSeq™ Analysis Report

